ACE AF-33

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Kort beskrivelse:

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Author liability:

All information herein was created conscientiously and under assumption of correctness. The creator cannot be held liable though for the use of any kind whatsoever and any ensuing consequences. The creator is not liable for any costs or expenses for damage that are the result of any information within this document. Any information is subject to change without prior notice.

Short description:

The ACE AF-33 is a digital electrochemical alcohol testing device, which measures the alcohol concentration in the breath and calculates the blood alcohol level from that. The device can be operated by only two buttons. The fuel cell sensor inside the ACE AF-33 can be calibrated. You can send the device to ACE Technik for a calibration service.

For each breath test a mouthpiece is needed. You can interpret the test result 1:1 (e.g. the display shows 0.36 = 0.36 per mill).

Scope of delivery:


- Digital electrochemical breathalyzer ACE AF-33
- 6 mouthpieces
- 2 pieces AAA batteries (1.5 V)
- Carrying case
- Operating manual

Introduction:


Thank you for purchasing our breathalyzer ACE AF-33. It measures alcohol concentration in human breath using an electrochemical sensor and offers reliable breath alcohol test results. The test results have got a tolerance of +/- 0.05 per mill. The device must only be used on basis of the instructions mentioned in this manual.

Label of individual components:

FRONT

BACK

This way you receive a result as accurately as possible:

- a) Wait at least 15 minutes after your last sip of alcohol. Do not eat, drink or smoke during that time.
- b) Do not open the casing.
- c) Do not keep the ACE AF-33 in cold, hot, moist or badly ventilated places.
- d) Keep the ACE AF-33 out of the reach of children.
- e) Conduct three inlet tests after long storage before conducting the actual alcohol test.
- f) Do not let the ACE AF-33 fall.
- g) Change batteries after battery alarm. Diminished battery capacity leads to incorrect test results.
- h) Do not conduct tests in loud (sound waves disrupt the sensor), draughty or badly ventilated places.
- i) After a test result of 1.00 per mille or higher the sensor needs a resting time. In this case wait at least three minutes before the next test.
- j) Do not conduct tests in very cold, very hot or unventilated, smoky ambience.
- k) Store the ACE AF-33 in a way to avoid penetration by dust or liquids.
- l) Blowing too hard may cause damage to sensor as saliva flows into sensor.

Preparation for the test:

1. Install 2 pieces of AAA batteries (alkaline, 1.5 V) in the battery compartment. Remove the battery compartment cover at the back of the device carefully and insert the 2 pieces of AAA batteries. Mind the battery polarity in the process (+/-).
2. For the active mode operation put a mouthpiece on the mouthpiece connection, which is on the side of the device. Make sure it is attached properly. Use a different mouthpiece for each test person. Avoid blowing cigarette smoke into the tester!


ACE Technik recommends:

- Always use a new mouthpiece for each person and each test which is conducted. ACE Technik assumes no liability for any damages which arise when several people use the same mouthpiece for multiple tests..
- Don't use the enclosed but the hygienic mouthpieces from the ACE series! During each test the ACE AF-33 internal and external gets in contact with human secretion (saliva, blood, sweat etc.)!
- Use safety gloves when pulling out a mouthpiece which was used by another person than yourself! You can buy appropriate gloves at ACE Technik.

Conducting a measurement:

1. Press power button. The total number of tests performed will be displayed.

2. The remaining recalibration date will be displayed after total number of tests. The Blue LED on both side of corner will flash with the message of "CAL" on the screen when the countdown reaches 0 day.

3. The countdown will be started automatically, and begin to inhale a deep breath when the countdown approaches 2.

4. When the countdown reaches zero with a double beep sound, three dashes will be displayed and „START“ will flash. Blow through the mouthpiece for five seconds until there is a double beep sound.

5. When the sampling is successful, LCD outer lines will be turning along the edge of LCD as you can see in the left picture.

6. The estimated BAC (Blood Alcohol Concentration) will be displayed after a few seconds (depending on the alcohol level) and the readings will flash for several seconds.

If you want to test again, press the power button to restart the countdown cycle.

- The tester will be turned off automatically after Step 6 in 15-20 seconds.
- The forced power off will be available in Step 3 and Step 4. Press and hold power button for 2-3 seconds to turn tester off.

Important notes:

When you have eaten, drunken alcohol or smoked a cigarette, wait at least 20 minutes before starting the measurement. If the test result is below 0.07 per mill then the device shows 0.00 per mill on the display. This happens because at a low range of alcohol concentration, there is a possibility to measure other material in your mouth.

Only the usage of original ACE mouthpieces

- guarantees the accuracy of the test results
- guarantees 100 % safety of hygiene
- guarantees full warranty of this device

ACE Technik warns:

Even if you get a result with the meaning „Safe range“, you should not drive after having consumed alcohol or do any other activities which are badly affected by the consumption of alcohol.

Memory mode:

The tester will recall up to 10 test results and their measurement time. To enter the memory mode, press memory button for about 2 seconds during the countdown (Step 3). The most recent test result will appear - first the measurement time will be shown (in three display methods, see below), then the test result will appear right after the measurement time. To see the next saved test result (and before its measurement time), just press the memory button sequentially.

The display method of measurement time:

1. Measured time up to 600 minutes ago will be shown as “-elapsed time”. For example:
Measured 1 minute ago -> -001
Measured 23 minutes ago -> -023
Measured 523 minutes ago -> -523
2. Measured time from 601 minutes ago to 1 day will be shown as “d000”. For example:
Measured between 601 minutes ago and 24 hours ago -> d000
3. Measured time from 1 day ago to 999 days ago will be shown as “d+ elapsed days”. For example:
Measured 1 days ago -> d001
Measured 23 days ago -> d023
Measured 523 days ago -> d523

The tester will recall up to 10 test results with measurement time. If the memory is full, the oldest test result including measurement time will be deleted automatically.

Notes to the battery law (BattG):

The device for which this operating manual is contains two AAA alkaline batteries (1.5 V) in delivery condition. Those batteries are not made for recharging.

Trying to recharge can cause damage to property or persons. Also the use of rechargeable batteries is not possible with this device, because it can cause damage to the device.

When the batteries are depleted, they must not be thrown in the household waste. Used batteries may contain harmful substances that can cause damage to the environment or health. Please return batteries / rechargeable batteries to the store or to the buyback centre. The return is free of charge und mandatory by law. Please throw only depleted batteries into the containers and tape over the poles of lithium batteries. For battery exchange only pull off the battery compartment cover and change the depleted batteries for new ones. Mind the polarity: The negatives poles of the batteries are in contact with the springs in the battery trough. Slide the battery compartment cover over the casing until it snaps in. All batteries and rechargeable batteries will be recycled. This way, precious resources like iron, zinc and nickel are gathered. Battery recycling is the easiest lesson in environmental.

The garbage bin means: Batteries and rechargeable batteries must not be thrown into the household waste. The signs below the garbage bins stand for:

- Hg: Battery contains mercury
- Cd: Battery contains cadmium
- Pb: Battery contains lead

Calibration:

To maintain the accuracy of the ACE AF-33, you can send the device for a calibration to ACE Technik. The calibration is not mandatory, but the accuracy of ACE AF-33 will decrease steadily without it. Send the device with the note 'Alkomat calibration service' to the company address stated on the invoice you received with the device. We recommend sending it as a parcel, so that in cases of loss an inquiry or proof of delivery is possible.

Due to changes in the product, the procedures stated can change. This can take place without prior notice or adaptation of the operating manual and is no cause for complaint, return or price reduction. The user is responsible to recognise deviations between the operating manual and the actual handling of the ACE AF-33 and he has to react accordingly.

Display

Meaning/Action

Full battery capacity.

Empty battery capacity: Warns when battery is empty. The unit can be operated but it is recommended to install new batteries for correct results.

The batteries are empty, no tests are possible, replace the existing batteries with new batteries immediately.

When the sampling volume or blow power is not enough to test. Press power button again to repeat the measurement.

If there is no breath sample blown within 30 seconds, the LCD will display „OUT“. Press power button to restart the countdown cycle and retest.

When the operating temperature shown in the specification is not within proper range. Repeat the test after the unit has reached a temperature within the operation range.

Display

If the unit displays 'Err2' after Step 5, it means the excessive saliva went into the sensor during the test. Wait 20 minutes and then try again.

When the total number of tests is over 500 times, the recalibration service is needed. Please send the device to ACE.

If the unit displays 'CAL' after countdown, it means the tester has been used over 365 days from its last calibration. Please send the device for a calibration service to ACE.

Remaining days to recalibration, e.g. d365 means 365 days left to recalibration. Contact ACE for proper recalibration service.

The elapsed day from the required recalibration date, which is a reminder to do the recalibration, e.g. d.001 means 1 day elapsed from the required recalibration date. Please send the device for a calibration service to ACE.

Technical data

Statement of the BAC	0.07 ~ 5.00 ‰ (fixed basic setting), under 0.07 ‰ the test result is shown as 0.00 ‰
Accuracy	+/- 0.05 ‰ at 0.5 ‰
Warm-up time	15 seconds up to max. 2 minutes
Response time	Within 12 seconds
Sensor	Fuel cell sensor
Energy supply	2 pieces AAA alkaline batteries (1.5 V)
Battery capacity	Approx. 1.000 tests
Dimensions of the device	47 mm x 108 mm x 17 mm (W x H x D)
Weight	85 g (including batteries)
Calibration	Every 6-12 months or after approx. 500 tests
Ambient conditions	Operation: 5 - 40° C Storage: 0 - 50° C

Warranty and liability

Warranty

The ACE Handels- und Entwicklungs GmbH will give a warranty of two years after purchase for material and processing errors for the product, to which this warranty and liability declaration is attached, in all versions and designations; a warranty of one year for the sensor. This warranty includes the electronics as well as the casing or the sensor without sampling components. Components that are inevitably destroyed upon effectiveness (such as fuses) are not part of this warranty. The warranty is valid from the date of delivery or from the date of first activation (always the shorter period), if this is not stated otherwise in the accompanying literature of the product by ACE Handels- und Entwicklungs GmbH.

Limitation of warranty and liability

1. Product description

The ACE Handels- und Entwicklungs GmbH rejects any further expressly stated or implied warranty, including the marketability and suitability for a determined purpose. This warranty especially excludes calibration. After changing the sensor, calibration is absolutely necessary.

When an exchange of the sensor or a calibration is necessary for a product, the user may receive an estimate of costs. Should the user refuse the estimate, the ACE Handels- und Entwicklungs GmbH will send the product back unrepaid / uncalibrated and the product will be labelled accordingly.

Should the product not fulfil the warranty, the claim of the buyer and the obligations of the ACE Handels- und Entwicklungs GmbH are limited exclusively and at the discretion of the ACE Handels- und Entwicklungs GmbH to an exchange or repair of the affected goods or to a reimbursement of the purchasing price for the faulty goods.

2. Application restriction

The ACE Handels- und Entwicklungs GmbH is on no account liable for additional damage, incidental damage or consequential damage such as, among others, lost profit and downtime that result from selling, producing or using products, which are part of this application restriction. This is independent of any claims,

such as from contracts, negligence or other culpable actions.

The product, which is accompanied by this warranty and liability declaration, no matter in which version or designation, can especially not be used legally binding to determine capability to operate vehicles or to execute action that could lead to damage to oneself or others under the influence of alcohol.

The thorough check for damage of all products upon receipt, a calibration suitable for the given use by the buyer and the proper use, maintenance and service according to the instructions of the accompanying product literature by the ACE Handels- und Entwicklungs GmbH are mandatory prerequisites for this warranty. Maintenance and service work by non-qualified persons, as well as using non-approved consumables or spare parts lead to an expiry of this warranty. As for all technically challenging products, it is a fundamental prerequisite and mandatory for the warranty by ACE Handels- und Entwicklungs GmbH, that all users of these products are familiar with the properties and restrictions of the product as stated in the respective literature.

3. Advice restriction

The buyer confirms that he has checked the intended use and suitability of purchased goods for any purpose in sole responsibility. Between the two parties it is expressly agreed upon that all technical or other advice by ACE Handels- und Entwicklungs GmbH regarding use of the goods or services is free of charge and was made with to the exclusive risk of the buyer. Therefore no obligation or liability can arise to the ACE Handels- und Entwicklungs GmbH from any recommendation or their result.

Please note:

The exact result of a person's blood alcohol concentration (BAC) can only be determined through a blood test. All indicated values determined through a breath test with this device are a clue and not legally binding. This device must be used strictly in accordance with the instructions given in the user manual, which is delivered with the device.

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ACE Handels und Entwicklungs GmbH